

The Engineering Division of YMC Europe GmbH  
 Schoettmannshof 19  
 D-46539 Dinslaken, Germany  
 phone: +49 (0) 2064 / 427-0, fax: +49 (0) 2064 / 427-222, E-mail: [info@kronlab.com](mailto:info@kronlab.com)

## Questionnaire for Kronlab Pilot Columns\*

*Please note: not all combinations are possible*

↻ **inner diameter (mm):** \_\_\_\_\_  
 ↻ **min. bed length (mm):** \_\_\_\_\_  
 ↻ **max. bed length (mm):** \_\_\_\_\_  
 ↻ **max. pressure (bar):** \_\_\_\_\_  
 ↻ **material column body:**      Glass       Acrylglass       SS   
 ↻ **O-ring material:**              EPDM       Viton®       Kalrez®   
 ↻ **biocompatible:**                  yes                   no   
 ↻ **frit material:**                    Polyethylene       SS       Titanium   
 ↻ **porosity (µm):**                    \_\_\_\_\_  
 ↻ **capillary connections:**              Polyethylene               PTFE               SS   
 ↻ **capillary outer diamter:**              1/8"                       1/4"                       1/2"   
 ↻ **capillary connection with valves:**                      yes                       no   
 ↻ **heating/cooling jacket:**              yes                       no   
 ↻ **temperature range:**                    \_\_\_\_\_  
 ↻ **column stand**                      Fixed stand       wheeled trolley       direct of plant   
 ↻ **bubble trap**                          yes                       no   
 ↻ **packing device:**                      yes                       no   
 ↻ **packing material:**                    \_\_\_\_\_  
 ↻ **solvents:**                              \_\_\_\_\_

Customer's address:	
name: _____	phone: _____
company/institute: _____	fax: _____
department: _____	email: _____
address: _____	
country: _____	